

Magnet Summary

Manufact.	Mode	Tone	RRM	BOL	ERI	Change	Program.
Medtronic	D/V/AOO	No	Off	85	65	Binary	Yes**
Biotronik	D/V/AOO	No	Off	90	80	Binary	Yes
Sorin/ELA	D/V/AOO	No	Off	96	80	Gradual	No
St Jude**	D/V/AOO	No	Off	98.6/100	86.3/85	Gradual	Yes
Boston Sci	D/V/AOO	No	Off	100	85	Gradual	Yes

Biotronik: Asynch--90 bpm continuously, Auto—90 for 10 beats only, Sync—no change

St Jude: Battery Test, Battery Test + EGM, OFF, OFF + EGM

Bost Scien: Pace ASYNC, EGM, OFF

Medtronic Pacers that will not respond to magnet for 60 minutes with standard end session: Kappa, EnPulse, Adapta, Versa, Sensia

St Jude Microny, Regency, Accent, Nuance, Anthem, Assurity, Endurity: 100/85